



Immunostaining and Super-Resolution Imaging of Structures In Cardiac Myocytes

> Shirley Zhang Auckland, New Zealand July 20, 2012

Proposed Research




 Super-resolution fluorescent imaging of ryanodine receptors (RyR) and microtubule structures in mouse cardiac myocytes with altered expression of protein junctophilin-II

 Compare relative positions of labeled structures to wild type samples

Progress This Week



 High resolution imaging of TAC control mouse cells for ryanodine receptors (RyR) and junctophilin-II (JPH)





(left) RyR in green Alexa680, JPH in red Alexa750; (right) RyR in red, JPH in green; Switching the secondary antibodies for the two structures seemed to show better results


Progress This Week



 Comparison of TAC1 and TAC2 strains of cells showed no significant differences







Future Plans



 Investigate effectiveness of antibody cocktail in staining tubule system in cells

- Sample 1: Na/Ca exchange only
- Sample 2: Caveolin 3 antibody only
- Sample 3: Both N/C ex and Cav 3

Cultural Aspect







Political protests against the New Zealand government, and street performers downtown

Cultural Aspect





My first bottle of L&P, a famous drink of New Zealand





Many thanks to:

- Gabriele Wienhausen, Teri Simas, Peter Arzberger UCSD PRIME
- Masahiko Hoshijima UCSD Mentor
- Christian Soeller Host Mentor
- The University of Auckland School of Medical Sciences
- The National Science Foundation



